



CENTRA Active

Made for active living.

www.siemens.com/hearing

SIEMENS



There's no reason a hearing loss has to slow you down. Now there's CENTRA Active™ from Siemens. It's water-resistant. Rechargeable. And offers the world's most advanced hearing technology – all in a new, comfortable, discreet design that's easy to wear. So you can do all the things you love to do without missing a beat. That's the beauty of CENTRA Active. For active people, it's central to what really matters.


Comfort, simplicity and style.

So you can focus on enjoying your life.

We designed CENTRA Active for easy operation and maximum comfort. It's completely automatic, so there are no controls to fuss with. And its unique receiver-in-the-canal design makes it tiny, easy to wear, and extremely comfortable.

CENTRA Active is ergonomically shaped for optimum placement in the ear canal and is virtually invisible when worn. Fitting is fast and easy with a wide selection of open tips to fit all ear canals. You can also choose from a broad range of colors to match hair and skin tones.





Our exclusive C-Guard barrier protects the dome tip against ear wax and humidity.

Water rolls right off the instrument's housing thanks to nanocoating technology.

Don't sweat it.

Water-resistant with AquaProtect, so you can do the things you love to do.

With CENTRA Active, you don't have to worry about working up a sweat playing golf, exercising, or taking a stroll on a warm afternoon. CENTRA Active's exclusive AquaProtect™ technology resists sweat, moisture, and humidity.


CENTRA Active's housing is nanocoated – a proprietary Siemens technology that repels water and prevents moisture and debris from leaking into the instrument. We also added a GORE™ clip-on microphone cover that protects against damage from perspiration, water, and other liquids. And our exclusive C-Guard™ barrier protects the dome tip against ear wax and humidity for months at a time.

Bye bye batteries.

Rechargeable, so you never have to hassle with batteries.

Batteries have long been a bother for hearing instrument wearers. You always have to make sure you have them on hand. If not, that means another trip to the store. And even when you do have them, they're not the easiest things for some people to handle.

But that's not a problem for CENTRA Active wearers. CENTRA Active's stylish, easy-to-use charger means no more worries about buying and handling batteries. Simply place the instruments in the charger – up to two instruments at a time. After just five hours of charging, the charger shuts itself off and the instruments are ready to go all day and evening. For even greater convenience, CENTRA Active will also run on regular hearing instrument batteries.



The charger selects the recharge time automatically – approximately 5 hours for a fully-drained rechargeable battery.



The optional ePocket remote control gives you easy, finger-tip control over virtually all instrument functions and read-out status.

Unrivaled hearing performance.

Only CENTRA Active has these breakthrough Siemens technologies.

CENTRA Active is more than rugged and rechargeable. Thanks to exclusive Siemens technologies, it gives you unparalleled hearing performance.

CENTRA Active's **SoundSmoothing™** feature reduces sounds like rustling paper and clanging dishes without affecting speech for better comfort and more natural sound. For binaural fittings, CENTRA Active's **e2e wireless™** synchronizes both instruments so you can locate sounds more precisely. This wireless feature also gives you easy, finger-tip control over virtually all instrument functions with the optional **ePocket™** remote. And CENTRA Active's unique **DataLearning™** technology learns the volume preferences you make with ePocket and automatically adjusts the instrument's volume for better performance.

Relax, it's automatic.

CENTRA Active surrounds you with more automated features than any other system for maximum ease of use.

Automatic Feedback Cancellation. Reduces or eliminates feedback or whistling that can occur with other hearing instruments.

Speech and Noise Management. Reduces background noise, even between syllables and words, and reduces noises like traffic and air conditioners for better hearing comfort and performance.

eWindScreen. With other hearing systems, wind can be bothersome as it blows across the instrument microphone. CENTRA Active's eWindScreen technology automatically adapts to windy conditions for maximum comfort.

Advanced Multi-Channel Directional Microphone System. This adaptive system automatically improves speech understanding when moving or competing sounds are present.

True 16-Channel Digital Signal Processing. Makes optimal use of your hearing and makes it easier for your Hearing Care Professional to precisely fine-tune the instruments to fit your unique needs.

Ready for anything.

Water-resistant, rechargeable, and offering unsurpassed hearing performance, CENTRA Active is ready for anything, so you can do all the things you love to do. It's the kind of innovation you'd expect from Siemens, a global technology leader whose 125 year commitment to advancing health care, energy, communications and transportation continues to transform the way we live.



The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice.

The required features should therefore be specified in each individual case at the time of conclusion of the respective contract.

© 06.2008 Siemens AG
Order No. A91SAT-00029-99C2-7600
Printed in Germany
104300/13715

Siemens Audiologische Technik GmbH
Gebbertstrasse 125
91058 Erlangen
Germany
Phone +49 9131 308 0

www.siemens.com/hearing